

ROTARY FILM EVAPORATOR

Utility Equipment for Rotary Film Evaporator

CIRCULATING CHILLER WATER BATH

- The Specifications of Chiller Water Bath is as below

Rotary Evaporator Model	Chiller MODEL	Temp. range	Cooling Capacity		KW	Tank Capacity	Process Pump	DRIVE MOTOR REQD H. P.
			Ton of Refrigeration	Kcal/Hour				
GRFE 5	FRAC 600	+5°C to +25°C	0.5	1500	1.7	20	35 LPM @ 0.3 Bar	0.25
GRFE 20 /GRFE 10	FRAC 1000	+5°C to +25°C	1	3000	3.5	40	50 LPM @1.2 Bar	0.5
GRFE 50	FRAC 2000 3P	+5°C to +25°C	2	6000	7.0	75	50 LPM @1.2 Bar	0.5
GRFE 100	FRAC 5000	+5°C to +25°C	5	15000	17.5	150	100 LPM	1

