

COUPLINGS

PTFE "O" RING WITH LOCKING COLLAR

These PTFE O rings are specially made to use as gaskets in glass fittings. These are provided with a collar which helps to locate it on the glass end correctly.

Cat.Ref.	DN	D	D1	d	L
TR0.5	12	26	18	3	5
TR0.7	15	28	17	3	5
TR1*	25	42	33	3	5
TR1.5*	40	57	48	3	5
TR2*	50	70	59	3	5
TR3*	80	100	88	3	5
TR4*	100	134	119	4	6
TR6*	150	186	168	4	6
TR8*	200	235	215	-	10
TR9*	225	260	236	4	7
TR12*	300	342	318	4	7
Above 12" TM used					

TM SIZE	OD	ID	THICKNESS
TM16	463	400	5
TM18	525	455	5
TM20	585	500	5
TM24	684	600	5
TM28	780	709	5
TM32	916	845	5

