




BOTTLES

2351

BOTTLES, ASPIRATOR WITH OUTLET FOR STOPPER

CAT. REF.	CAP. ML.	APPROX O.D. MM	APPROX HEIGHT MM	APPROX NECK I.D. MM	APPROX OUTLET I.D. MM	QTY PER CASE	Group
2351/32	5000	190	320	44	28	1	D
2351/35	10000	220	420	54	28	1	D
2351/40	20000	300	535	54	28	1	D



2356

BOTTLES, ASPIRATOR , WITH INTERCHANGEABLE STOPPER AND STOPCOCK

CAT. REF.	CAP. ML.	APPROX O.D. MM	APPROX HEIGHT MM	APPROX NECK I.D. MM	APPROX OUTLET I.D. MM	Qty Per Case	Group
2356/32	5000	190	320	45/40	29/32	1	D
2356/35	10000	220	420	55/44	29/32	1	D
2356/40	20000	300	535	55/44	29/32	1	D

